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VOCATIONAL TRAINING CORPORATION (VTC)

**TECHNICAL OPERATIONS IN THE PHARMACEUTICAL
PRODUCTION UNITS
(LEVEL 3)**

A COMPETENCY BASED APPROACH

**Equipment and
Production Rooms:
Cleaning and
Maintenance (C&M)**

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Acknowledgment

This publication has been made possible with the invaluable contribution by many stakeholders including competitive inputs and leadership by the pharmaceutical industry in Jordan, instructors, and consultants.

The Ministry of Labour and ETVET Council through Vocational Training Corporation (VTC) and the European Union with technical assistance from GOPA Worldwide Consultants have partnered with the following industry experts, the Jordanian Association for Pharmaceutical Manufacturers (JAPM) to develop the curriculum framework and the resulting learning materials.

The VTC would like to thank particularly the following Consultants and Members of the Curriculum Development Panel for their works in research, developing the content, review and finalising these materials.

Dr. Hanan Sboul	JAPM
Mr. Rami Dari	Ram Pharma
Eng. Sana' Abbas	JPM
Dr. Ayyad Rumman	Ayyad Rumman Pharma Consulting Centre
Mrs. Suzanne Grigoleit	VTC
Dr. Judy Saleh	PCOE - VTC
Eng. Hussain Zuhairi	PCOE - VTC
George Kolath	Lead ISTE, GOPA

No part of this publication may be reproduced, stored in a retrieval system, or transmitted in any form or by any means, electronic, mechanical photocopy and recording or otherwise, without the written permission of Vocational Training Corporation.

The publishers have made every effort, where required, to acknowledge copyright holders, but if any have been overlooked, the necessary arrangements will be made at the first opportunity. We welcome any information that will enable us to correct, in subsequent edition, any inaccurate or omitted credit.

No liability is assumed with respect to the information contained in.

Published by

Vocational Training Corporation
Amman - Jordan

With financial and technical assistance by the European Union and GOPA Worldwide Consultants respectively.

© Vocational Training Corporation
2014

All rights reserved.

This Project is Funded by the European Union

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Contents

- 10** **Overview of the Module**
- 14** **cGMP Requirements for Cleaning in Pharmaceutical Production**
- 22** **Prepare Cleaning Tools, Equipment and Materials**
- 34** **Cleaning of Equipment**
- 42** **Cleaning of Production Rooms**
- 50** **Documentation of the Cleaning Process**

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Overview of the Module


In the Modules on cGMP (P101) and OEHS (P102), we learnt about various aspects of quality management and quality assurance including safety requirements. While discussing quality management, we learnt that hygiene and cleanliness are very critical in the pharmaceutical industry in order to maintain low cost sustainable, high quality and safe operations in a manufacturing facility.

This module, Equipment and Production Rooms: Cleaning and maintenance focuses on the critical importance of maintenance and sustainable production procedures to be followed by the pharmaceutical industry as per cGMP norms and

international regulations and generally accepted industrial standards.

This Module is essential for your employability as a Technician in the Manufacturing Industry. At the end of this module, you will be able to demonstrate knowledge, skills and competences in the SOPs for cleanliness and simple maintenance of pharmaceutical production rooms and equipment.

These Units, when completed, will provide you with the necessary competencies to carry out and document effective cleaning of the equipment and production rooms in the pharmaceutical industry.

Learning outcomes

When you complete the Module, you will:

- Illustrate procedures for cleaning pharmaceutical manufacturing equipment, clean rooms and machineries.
Describe and correctly use various detergents and disinfectants in the cleaning process.
- Apply safety measures to be included while using cleaning tools, equipment and materials.
- Demonstrate compliance with recording and reporting procedures.

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Credit value

The credit value of this module is **3 credits** which equals to **65 hours** of nominal learning time. This includes theory lesson (face-to-face), supervised and self-directed practice sessions, and continuous and summative assessment (learning controls).

Assessment

At the end of every Unit there are learning and self-assessment activities. The learning and assessment activities provide you with the opportunity to apply knowledge and skills you have acquired. You should successfully complete all the learning and assessment activities as evidence that you actually mastered the skills and related knowledge in this Module.

Your Instructor too will carry out a formal assessment as evidence of your certification. As before, the mode of assessment will be both continuous (formative assessment) and summative (external assessment).

In this Module, your continuous assessment, both theoretical and practical consisting of 30% of all Units, will be considered for the final grading.

The final assessment of the entire Module will be an examination of both theory (20%) and demonstration or practicals (80%), which will be carried forward for the purpose of certification.

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Unit 1

cGMP Requirements for Cleaning in Pharmaceutical Production

Learning outcomes

At the end of this unit, you will:

1. Understand the concept of cleaning in pharmaceutical industry
2. Interpret the requirements of cGMP in cleaning equipment, cleanrooms and clean areas
3. Apply your skills to comply with the regulatory requirements

1.1 Why to clean?

Clean and healthy is a saying you may have come across at your school days. Air pollution and other contaminations make us sick. Imagine then how important it is to have safe medicine to treat sicknesses or diseases!

The cleaning of equipment and production is an important aspect of pharmaceutical production. They are very critical as it is a regulatory issue in as much as afford remarkable savings due to avoidance of breakdowns or lower batch production failures.

The cGMP (WHO 12.1 and 12.2; USFDA 211.182) requires that **equipment, production areas and machineries must be cleaned regularly on set standards to avoid cross contamination and degeneration of the drug.**

Through the Units in this module, your Instructor will help you through a number of theoretical and

practical sessions on various cleaning steps, appropriate disinfectants used in the cleaning process, and also their accurate and timely documentation.

Definitions and meaning

Cleaning

Removal of all remnants from surfaces, floors, machines and machine parts, equipment and devices.

Sterilisation

Extermination of all germs (micro-organisms) from surfaces and materials. This can be done using either the physical methods (heat, ration, percolation, or ultrasonic waves), or chemical methods (sterilising solutions or gases).

Disinfection

Extermination of the vegetative state of bacteria rather than the endospores. Disinfection kills viruses or at least decreases their viral activity. Antiseptics are used for living tissues, and disinfectants are used for non-living objects.

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1.2 Importance of effective cleaning and sanitization process

A formal cleaning and sanitisation process is essential in the manufacturing of pharmaceutical products, and specific and written internal procedures, as we have pointed out, should be strictly followed for cleaning and sanitation, because it:

- Assures microbiological quality ('purity') of the products.
- Meets legal regulations were required.
- Minimizes the chances of degeneration of products contributed to processing, packaging and storage.
- Avoids the cost associated with failure of products.
- Helps maintain the organization's commitment to quality.
- Manintain good reputation of your company.

Factors that affect efficiency of sanitizers of disinfectants may include the following:

1. Time- the duration of exposure to sterilizer or disinfectant.
2. Temperature.
3. Concentration, method and date of preparation.
4. The microorganism.
5. The material to be sterilized.
6. Level of cleanliness of the material to be sterilized.
7. Cleaning methods used.

1.3 cGMP requirements and SOPs for the cleaning of equipment and production rooms

cGMP guidelines demand that the facilities, the buildings and surroundings need to be accessible and easy to clean. In stricter terms, they call for stringent environmental controls, beginning with planning and construction of rooms. In order to achieve clean air standards rooms must be separate, adequately-sized and have equipment to control humidity, dust, air pressure, temperature and microorganisms. (we will later explore the environmental monitoring programme.)

What are cleanrooms?

Pharmaceutical facilities are made- up of a series of rooms called cleanrooms.

Cleanroom is controlled placement where different products are manufactured.

In the pharmaceutical production, the very survival of the manufacturer depends on the safety of the finished products. The basic function of a cleanroom, therefore, is to protect the manufactured product from contamination.

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The importance of clean standards is also stressed by the cGMP stating that “any building used in the manufacture, processing, packing or holding of a drug product shall be maintained in a clean and sanitary condition.”

You will have realised that every pharmaceutical manufacturing plant should have well laid out SOPs for sanitation which contain the following area-wise levels of cleanliness and sanitation:

- Aseptic (germ free/sterile) area
- Processing areas
- Packing areas
- Storage areas and utility areas
- Change rooms
- Offices and corridors
- Usage procedures of detergents and disinfectants

- Cleaning techniques in detail.
- Contact times.
- Rinsing.
- Frequency of cleaning and disinfecting.
- Procedure for transfer of disinfectants and cleaning agents in and out of the clean areas (including the procedure for sterilisation of disinfectants).
- Holding time for detergents and disinfectants.


Fig 1.1 Diagram of a clean production facility (Source: Tuttnauer)

1.3.1 The SOPs on cleaning the production rooms and equipment should describe:

- The type of detergents and disinfectants to be used (which are compatible for different equipment and processes).
- The frequency of rotation of disinfectants.
- A list of suitable cleaning materials.

1.4 Environmental quality monitoring


The effectiveness of the cleaning and sanitisation process should be accessed through an environmental monitoring programme.

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Environmental monitoring is a programme that examines the numbers and occurrence of viable micro-organism and non-visible particles like dust, dry skin cells etc.

Ideally, environmental monitoring is targeting those areas where the risk of contamination cannot be controlled easily.

Individual assignment

Fill in appropriate details in the following report format and review them with your instructor.

JORDAN PHARMACEUTICALS Pvt Ltd.

AREA CLEANING RECORD- RAW MATERIALS, PREVENTIVE MAINTENANCE & FINISHED GOODS

Department:						Month:					
Area:						Ref. SOP No:					
Room No:											
Date	Floors	Furniture	Office Equipments	Containers	Pallets	Racks	Standards wt	Cleaned by:	Time of cleaning		Signed By Officers
									From	To	

Table 2.1 Area Cleaning Record Format (Daily)

Learning activities

1. Assisted by your Instructor, design an SOP on how to disinfect a production room. Further to this, undertake a field visit to a pharmaceutical firm, with prior appointment through your Instructor, to interview the Production Manager on various cleaning processes at a production room. Upon your return, review the actual working procedures versus your initial SOP.
2. Carry out further reading on disinfectants and cleaning agents in the pharmaceutical industry.

Test yourself

1. What is cleaning: at the household, industry, pharmaceutical industry level?
2. Explain in brief the requirements of the cleaning SOPs in pharmaceutical industry.
3. Design a SOP for cleaning the Cleanroom.

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Explore: <http://www.ilo.org/safework/countries/arab-states/jordan/lang--en/index.htm>
en.wikipedia.org/wiki/Occupational_safety_and_health

Unit 2

Prepare Cleaning Tools, Equipment and Materials

Learning outcomes

At the end of this module, you will:

1. Differentiate between various cleaning materials and tools.
2. Understand various procedures and hazards in dealing with cleaning materials.
3. Prepare cleaning materials and equipment.

2.1 Introduction

In Unit One, we learnt that clean rooms and production areas are to be regularly cleaned and disinfected. In this Unit, we focus on various types of cleanings, substances, and tools used to clean production rooms and equipment.

cGMP guidelines demand that the cleaning and sanitisation programme in the pharmaceutical industry should recognise specified cleanliness standards and control microbial (contagious) contamination of products. It should also aim at preventing the chemical contamination of pharmaceutical ingredients, product-contact surfaces and/or equipment, packaging materials, and ultimately the drug products.

You will also learn to use, during your practical sessions, the most common cleaning equipment used for physical removal of dust and dirt such as vacuum cleaners, brooms, brushes, floor mops, sponges and jet cleaners.

Let us now briefly understand the various definitions of terminologies that you should be very familiar as they are closely related to the cleaning and sanitation programme in the pharmaceutical industry.

2.2 Definitions of terms - disinfectants and antiseptics

Antiseptic

An agent that inhibits or destroys microorganisms on living tissue including skin, oral cavities, and open wounds.

Chemical disinfectant

A chemical agent used on inanimate surfaces and objects to destroy infectious fungi, viruses, and bacteria, but not necessarily their spores.

Cleaning agent

An agent for the removal from facility and equipment surfaces of product residues that may inactivate sanitizing agents or harbor microorganisms.

Decontamination

The removal of microorganisms by disinfection or sterilization.

Disinfectant

A chemical or physical agent that destroys or removes vegetative forms of harmful microorganisms when applied to a surface.

Sanitising agent

An agent for reducing, on inanimate surfaces, the number of all forms of microbial life including fungi, viruses, and bacteria.

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Sporicidal agent

An agent that destroys bacterial and fungal spores when used in sufficient concentration for a specified contact time. It is expected to kill all vegetative microorganisms.

Sterilant

An agent that destroys all forms of microbial life including fungi, viruses, and all forms of bacteria and their spores. Sterilants are liquid or vapor-phase agents.

2.3 Types of disinfectants

A disinfectant is one of a diverse range of chemicals which reduces the number of micro-organisms present in an object.

Disinfectants are usually available in concentrated solutions. And you should remember that they are to be diluted as recommended by manufacturers and documented suitably.

Micro-organisms can be classified as: Bacterial spores, Mycobacteria, Nonlipid-coated viruses, fungal spores and vegetative molds and yeast, vegetative bacteria, and lipid-coated viruses.

There are many different types of disinfectants used in the pharmaceutical industry, with different ranges of activity and modes of action.

Disinfectants are selected for use based on mode of action, efficacy, compatibility, cost and current health and safety standards.

Detergent is a water-soluble cleansing agent which combines with impurities and dirt to make them more soluble, and differs from soap in not forming a scum with the salts in hard water.

Surface disinfectants are usually classified into the following categories:

- Non-oxidizing Disinfectants such as alcohols, amphoteric and penicillins.
- Oxidizing Disinfectants such as peracetic acid and hydrogen peroxide.

2.4 Types of cleaning equipment and tools

a. Mops- Mops may be made of sponge and should preferably be of 'squeeze' type. The mops should be cleaned in soap water after every use, rinsed from disinfectant or cleaning agent allowed to dry in air and put away in an airy place.

b. Clothes and sponges- They are very useful for working on surfaces, walls, ceiling and even on equipment, These are of absorbent materials and remember to store all clothes and sponges used in cleaning in a dry condition.

c. Vacuum cleaners- These are more effective to clear dust than brooms or brushes. Care must be ensured to fit them with suitable filters. Remember to use dedicated vacuum cleaners in the sterile areas. The dust collecting bag of the cleaner should be cleaned after every use (in the case of portable vacuum cleaner).

Disinfectants are antimicrobial agents that are applied to non-living objects to destroy microorganisms that are living on the objects.

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d. Mechanical floor cleaners

These are suitable for large manufacturing areas. They are very effective to remove sticky stains on walls and floors.

e. Jet cleaners- These devices spray hot water or detergent solutions, fitted with various types of nozzles to enforce the sprays.

f. Brooms- They are not generally recommended as brooms leave particles of dust on the floor or surfaces.

In groups, draw up a list of common disinfectants used in the cleaning/sterilisation processes within the industry. You should match this list with actual products/samples.

2.5 Types of cleaning

In the industry, we employ three types of cleaning methods to ensure highest levels of cleanliness and decontamination. They are **Full Clean, Wet Clean** and **Dry Clean Methods**. Let us briefly discuss them in detail.

I. Full clean method

It is a full removal of products remnants from production rooms, equipment, or machines. In this process, disinfectants and cleansing agents are used. This type is used in the following cases:

- Between 2 different products.
- New equipment, devices, and machines.
- After Fogging.

II Wet clean method

It is a removal of products remnants from production rooms, equipment, or machines using water and disinfectants. It is used in the following cases:

Between 2 products of similar composition, yet different concentrations (from more to less concentrated).

- Between 2 products similar in composition, yet different in colors.
- After long duration of work, in line with SOP.
- After a long period of time after a full cleaning procedure for rooms, equipment, or machines, in line with SOP.

III Dry clean method

It is the process of removal of products remnants from production rooms, equipment, or machines without using water or disinfectants. It is used in the following cases:

- Between 2 products of similar composition and concentration.
- Between 2 products of similar composition, yet different concentrations (from less to more concentrated).

2.6 Preparation of tools and cleaning products

- Tools to be used: Powder vacuum, plastic brush, towel, whipper, hose, and a pump for cleaning products.
- Materials to be used: Water, hot water, distilled water, 70% alcohol, isopropyl alcohol, disinfectant, and cleaning agents, as specified in the SOPs.

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2.6.1 Procedures to be followed when cleaning

When you begin the process of cleaning, firstly, ensure the following PPEs:

1. Put on protective medical gloves when preparing cleaning solution and cleaning.
2. Wear appropriate facemask.
3. Wear protective eyeglasses.
4. Disconnect all machines from electricity when cleaning them.
5. Wear cleaning uniform.
6. Wear cleaning shoe.

2.6.2 Cleaning methods

In all cleaning methods you are expected to follow strictly the cGMP documentation: specific SOPs for cleaning, and general safety procedures.

Now let us briefly discuss the three categories of cleaning methods employed at the pharmaceutical industry.

In **Full clean method**, you are asked to follow the given procedures below:

1. Remove all things related to previous products.
2. Disconnect electricity from machines and devices to be cleaned.
3. Bring a copy of SOP of cleaning.
4. Fill in Forms and Records for this procedure.
5. Remove all things related to previous batch (label, identification, powder, tablets, strips etc.).
6. Dis-assemble parts that need to be cleaned.
7. Vacuum and isolate remnants of previous product.
8. Prepare concentration of cleaning solution in line with SOP.

9. Clean machine, production rooms, equipment, devices, and disassembled parts using cleaning solution.
10. After using cleaning solution, wash with regular water, and then with distilled water.
11. Wash/ wipe disassembled parts, machines, and tools with 70% alcohol.
12. Dry using compressed air or a towel.
13. Wash the floor, walls, ceiling, and lights.
14. Disinfect the floor, and walls using appropriate disinfectant in line with SOP.
15. Put cleaning tools and materials back to its place in washing room after using it.
16. Ensure cleaning by direct supervisor.
17. Complete forms and records documentation.
18. Call in supervisor to ensure completion of cleaning procedure.

In the **Wet clean method**, you follow all steps of full cleaning, except using cleaning solution in steps (8) and (9).

In the **Dry clean method**, you are asked to follow the given steps below:

1. Remove all things related to previous batch.
2. Disconnect electricity from machines and devices to be cleaned.
3. Bring a copy of SOP of cleaning.
4. Fill in Forms and Records for this procedure.
5. Remove all things related to previous batch (label, identification, powder, tablets, strips ...etc.).
6. Dis-assemble parts that need to be cleaned (if required).
7. Vacuum and isolate remnants of previous product.

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8. Clean the floor in a dry method.
9. Ensure cleaning by direct supervisor.
10. Complete forms and records of documentation.
11. Call in supervisor to ensure completion of cleaning procedure.

Remember

1. All records and documents should be filled from the beginning of work till the end of it.
2. Direct supervisor should be informed about any malfunctions that may occur during the cleaning process.
3. Cleaning wastes should be isolated, packed in appropriate canister, identified, and handed over to the official in charge, after it is being sorted into solid and packaging.

2.7 Equipment sanitisation

Effective cleaning and sanitization of equipment are important because equipment is not responsive to visual inspection and there is a good possibility of formation of biofilm on it.

The main method for cleaning industrial equipment is by making the mechanism for cleaning corresponding to the equipment itself. This can be achieved by use of pressure, heat, steam sterilization, mechanical removal, or chemicals, and is termed Clean-in-Place (CIP) or Steam-in-Place (SIP).

Prior to chemical or heat treatment, attempts must be made to remove process residues and particles using steam or high-pressure water cleaning.

Alkali-based disinfectants and detergents are commonly used for CIP systems, with sodium hydroxide among the most widely used. Such caustic alkalis can readily remove organic deposits without affecting the equipment.

2.8 Hazards and procedures for handling and use of different cleaning materials

cGMP 21 CFR 211.67, Equipment Cleaning and Maintenance, details the requirements for written procedures for cleaning, maintenance, and sanitization of pharmaceutical manufacturing equipment:

- These procedures should address the assignment of responsibility, establishment of schedules, details of cleaning operations, protection of clean equipment prior to use, inspection for cleanliness immediately prior to use, and maintenance of cleaning and sanitization records.
- Staff involved in disinfection requires training in microbiology, industry practices for cleaning and sanitization, safe handling of concentrated disinfectants, the preparation and disposal of disinfectants, and appropriate application methods.
- It should be emphasized that the preparation of the correct dilutions is critical because many disinfectant failures can be attributed to use of disinfectant solutions that are too dilute.
- Typically disinfectants used in aseptic processing and filling areas are diluted with Sterile Purified Water, and are prepared aseptically.

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- Alternately, the disinfectant may be diluted with Purified Water, and then sterile filtered to eliminate microorganisms that may potentially persist in a disinfectant. Diluted disinfectants must have an assigned expiration dating justified by effectiveness studies.

Because it is theoretically possible that the selective pressure of the continuous use of a single disinfectant could result in the presence of disinfectant-resistant microorganisms in a manufacturing area, some specialists have

Sporicidal Agents are substances used for killing spores.

suggested the rotation of disinfectants. It is prudent to augment the daily use of a bactericidal disinfectant with weekly (or monthly) use of a sporicidal agent. The daily application of sporicidal agents is not generally favoured because of their tendency to corrode equipment and because of the potential safety issues with chronic operator exposure. Other disinfection rotation schemes may be supported on the basis of a review of the historical environmental monitoring data. Disinfectants applied on potential product contact surfaces are typically removed with 70% alcohol wipes. The removal of residual disinfectants should be monitored for effectiveness as a precaution against the possibility of product contamination.

The greatest safety concerns are in the handling of concentrated disinfectants and the mixing of incompatible disinfectants. For example, concentrated sodium hypochlorite solutions (at a concentration of more than 5%) are strong oxidants and will decompose on heating, on contact with

acids, and under the influence of light, producing toxic and corrosive gases including chlorine. In contrast, dilute solutions (at a concentration of less than 0.5%) are not considered as hazardous.

Under no circumstances should disinfectants of different concentrations be mixed.

Material Safety Data Sheets (MSDS) for all the disinfectants used in a manufacturing area should be available to personnel handling these agents. Appropriate safety equipment such as face shields, safety glasses, gloves, and uniforms must be issued to personnel handling the disinfectant preparation, and personnel must be trained in the proper use of this equipment. Safety showers and eye wash stations must be situated in the work area where disinfectant solutions are prepared.

Learning activities

1. Identify, through internet search or books in your library, various cleaning and sanitisation materials used in the pharmaceutical industry.
2. Demonstrate and be conversant with application of cleaning materials

Test yourself

1. Classify different types of disinfectants used in pharmaceutical industry.
2. What are the hazards and precautions that need to be taken while implementing a cleaning process?
3. Why should disinfectants of different concentrations NOT to be mixed?

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Learning outcomes

In this Unit, you will:

1. Understand step-by-step procedure of cleaning production machines
2. Apply cGMP guidelines for cleaning production machines

Case study: Let us begin with a real case study serving as an example on the need for quality cleaning of equipment.

Cleaning of machineries forms a vital part of the production process

One event which increased US Food and Drug Administration (FDA) awareness of the potential for cross contamination due to inadequate procedures was the 1988 recall of a finished drug product, Cholestyramine Resin USP. The bulk pharmaceutical chemical used to produce the product had become contaminated with low levels of intermediates and degradants from the production of agricultural pesticides. The cross-contamination in that case is believed to have been due to the reuse of recovered solvents. The recovered solvents had been contaminated because of a lack of control over the reuse of solvent drums. Drums that had been used to store recovered solvents from a pesticide production process were later used to store recovered solvents used for the resin manufacturing process. The firm did not have adequate controls over these solvent drums, did not do adequate testing of drummed solvents, and did not have validated cleaning procedures for the drums.

Some shipments of this pesticide contaminated bulk pharmaceutical were supplied to a second facility at a different location for finishing. This resulted in the contamination of the bags used in that facility's fluid bed dryers with pesticide contamination. This in turn led to cross contamination of lots produced at that site, a site where no pesticides were normally produced.

Now in groups, as guided by your Instructor, you should:

1. Analyse the above case study; and,
2. Propose a step-by-step procedure how this contamination would have been prevented (you may recall here various types of cleaning you learnt in Unit 2).

3.1 Clening procedures

In Module 2, we discussed considerably on the necessity of clean environment and clean operations as critical to the pharmaceutical manufacturing of quality and effective drugs.

Cleaning procedure involves addressing such issues as- What products are being processed? What equipment and areas are involved? What are the materials of construction of the equipment and surfaces in the facility that require cleaning? What cleaning agents are being used in what concentration? and, What degree of cleanliness is expected by the regulations?

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In order to ensure a reliability and efficiency of cleaning, cleaning procedures should be sufficiently detailed (SOPs) to remove the possibility of any inconsistencies during the cleaning process.

You may track the following **parameters** during **cleaning procedures**:

A. Equipment Parameters to be evaluated

1. Identification of the equipment to be cleaned
2. 'Difficult to clean' areas
3. Property of materials
4. Ease of disassembly
5. Mobility

B. Residues to be cleaned

1. Cleaning limits
2. Solubility of the residues
3. Length of campaigns

C. Cleaning agent parameters to be evaluated

1. Preferable materials that are normally used in the process
2. Detergents available (as a general guide, minimal use of detergents recommended unless absolutely required)
3. Solubility properties
4. Environmental considerations
5. Health and safety considerations

D. Cleaning techniques to be evaluated

1. Manual cleaning
2. CIP (Clean-in-place)
3. COP (Clean-out-of-place)
4. Semi-automatic procedures
5. Automatic procedures
6. Time considerations
7. Number of cleaning cycles

3.2 The procedure of cleaning machines

Repeatedly we discussed in this Module on the need for cleaning equipment, utensils and production rooms. Primarily, this is carried out in order to prevent malfunctions and contamination.

There certain procedures for equipment cleaning and maintenance. Some of them, which you will come across at the manufacturing facilities, are:

1. Assignment of responsibility for cleaning and maintaining equipment
2. Maintenance and cleaning schedules, including sanitising schedules
3. A description of sufficient details of the methods, equipment and materials used in cleaning and maintenance of operations, and the methods of disassembling and reassembling equipment as necessary to assure proper cleaning and maintenance.
4. Removal of previous batch identification
5. Protection of clean equipment from contamination prior to use
6. Inspection of equipment for cleanliness immediately before use.

You will be given demonstrations and practical sessions on cleaning and sanitising production rooms and equipment with specific SOPs in this regard. Also, remember that you will be actually carrying out these operations (cleaning/sanitation) during specific modules such as dispensing raw materials, granulation/drying, and tableting or encapsulation.

The first step in the process of cleaning equipment is to identify the cleaning requirement of the machine. A piece of equipment or system will have one process for cleaning, however this will depend on the products being produced and whether the clean-up occurs between batches of the same product (as in a large campaign) or

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between batches of different products. When the cleaning process is used only between batches of the same product (or different lots of the same intermediate in a bulk process) the firm need only meet a criteria of, "visibly clean" for the equipment. Such between- batch- cleaning processes do not require validation.

3.3 Equipment design

Ideally, equipment and machineries used in the production process of pharmaceuticals should have a Clean-in-Place (CIP) system in them. When such a system is in place, the operator of the system should be adequately trained to clean such systems. The cleaning process for such machines should be in writing and validated regularly. In larger system where employees form large transfer lines, the flowchart of the assembly system must be analysed to identify critical checking points in the system and explicit written instructions should be provided on cleaning and validating the critical points.

Clean-in-place (CIP) is a method of cleaning the interior surfaces of pipes, vessels, process equipment, filters and associated fittings, without disassembly.

Whether or not CIP systems are used for cleaning of processing equipment, microbiological aspects of equipment cleaning should be considered. This consists largely of preventive measures rather than removal of contamination once it has occurred. There should be some evidence that routine cleaning and storage of equipment does not allow microbial proliferation.

Subsequent to the cleaning process, equipment may be subjected to sterilization or sanitization procedures where such equipment is used for sterile processing, or for non-sterile processing where the products may support microbial growth. While such sterilization or sanitization procedures are beyond the scope of these learning materials, it is important

to note that control of the bio-burden through adequate cleaning and storage of equipment is important to ensure that subsequent sterilization or sanitization procedures achieve the necessary assurance of sterility.

3.3.1 Sampling of cleaning

There are two general types of sampling that have been found acceptable. The most desirable is the direct method of sampling the surface of the equipment. Another method is the use of rinse solutions.

What is sampling?

Sampling may be defined as the process of removing appropriate number of items from a population in order to make inferences to the entire population.

- a. **Direct surface sampling** - Determine the type of sampling material used and its impact on the test data since the sampling material may interfere with the test. For example, the adhesive used in swabs has been found to interfere with the analysis of samples. Therefore, early in the validation program, it is important to assure that the sampling medium and solvent (used for extraction from the medium) are satisfactory and can be readily used.

Advantages of direct sampling are that areas hardest to clean and which are reasonably accessible can be evaluated, leading to establishing a level of contamination or residue per given surface area. Additionally, residues that are "dried out" or are insoluble can be sampled by physical removal.

- b. **Rinse samples** - Two advantages of using rinse samples are that a larger surface area may be sampled and inaccessible systems or ones that cannot be routinely disassembled can be sampled and evaluated.

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A disadvantage of rinse samples is that the residue or contaminant may not be soluble or may be physically sealed in the equipment.

Limits should be set for acceptance level of residue on equipment. Although there is no basic standards set by any authority because of the wide variety of equipment being used for bulk pharmaceutical production and finished goods, the firm should have a rationale for such acceptance standards.

Learning activities

Demonstrate use of cleaning machines in real life pharmaceutical setting.

Test yourself

What are the major factors that need to be considered while cleaning equipment in pharmaceutical firms?

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Learning outcomes

You will:

1. Demonstrate procedures of cleaning clean rooms and production areas
2. Apply cGMP guidelines for cleaning clean rooms and production areas

4.1 Introduction

Pharmaceutical manufacturing takes place within a series of specially controlled environments, known as cleanrooms. A cleanroom is a room in which the concentration of airborne particles is controlled to specified limits.

On the one hand, a cleanroom or clean zone for production is simply an area that is clean in terms of both particle counts and microbial counts. In addition, regulatory requirements for cleanrooms are detailed by EU GMP or the FDA guidelines.

The purpose of cleaning and disinfection is to achieve appropriate microbiological cleanliness levels for cleanroom for a suitable period of time. One important step towards achieving microbial control within a cleanroom is the use of defined cleaning techniques, together with the application of detergents and disinfectants.

The detergents and disinfectants used in pharmaceutical grade cleanrooms need to be of a high quality and effective at killing micro-organisms. Both correct product selection and cleaning

Sources of contamination

This is a partial list of some of the commonly known contaminants that can cause problems in some production room environments. Many of these contaminants are generated from five basic sources. The **facilities, people, tools, fluids** and the **product being manufactured** can all contribute to contamination. Review this list to gain a better understanding of where contamination originates.

Facilities

Walls, floors and ceilings; paint and coatings; construction materials; air conditioning debris; room air and vapors; and, spills and leaks.

People

Skin flakes and oil; cosmetics and perfume; spittle; clothing debris; and, hair.

Tool or equipment-generated

Friction and wear particles; lubricants and emissions; vibrations; and, brooms, mops and duster

Fluids

Particulates floating in air; bacteria, organics and moisture; floor finishes or coatings; and cleaning chemicals; and, deionized water

Products generated

Packaging debris; powder/granule particulates; and leftovers of tablets or capsules etc.

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techniques are important, particularly in relation to some of the newer cleanroom technologies.

4.2 Cleaning agent selection

The types of detergents and disinfectants used represent an important decision for the pharmaceutical manufacturer. There are various different types of disinfectant with different spectrums of activity and modes of action.

Pharmaceutical manufacturers are expected to use at least two disinfectants of different modes of activity in order to conform to current Good Manufacturing Practices. The use of a sporacidal disinfectant is recommended for sterile areas on an occasional basis, even where such a disinfectant does not form part of the standard set. When selecting disinfectants, it is prudent to opt for manufacturers who offer a range of disinfectants of different modes of activity (which target and kill microbial cells in different ways).

For a disinfectant to work effectively, "soil" (such as grease and dust particles) must be removed first using a suitable cleanroom grade detergent. Detergents are typically either sterile, neutral solutions or have added cleaning chemicals, such as premium sterile non-ionic surfactants.

When selecting a detergent a check should be made to determine that it is compatible with the disinfectant because some detergents, where there are residues, can neutralise the active ingredient in certain disinfectants. This can be overcome where detergents are purchased from the same manufacturer as the disinfectant.

4.3 Cleaning techniques for cleanroom

The cleaning technique used for cleanrooms should be defined and standardised. In most cases, cleaning of Cleanroom is applied with the three types of cleaning we have already discussed in Unit

2.4 (Full Clean, Wet Clean and Dry Clean Methods). However, in this subsection, let us review in brief, some of the cleaning techniques to be employed by the Cleanroom Operators.

General cleaning of the cleanroom can be demonstrated through the following steps:

- **Sweeping away dust and debris (if applicable).**
- **Applying a detergent solution through wiping or mopping.**
- **Ensuring that the detergent has dried.**
- **Applying a disinfectant solution through wiping or mopping.**
- **Keeping the surface wet until the contact time has elapsed.**
- **Removing disinfectant residue through wiping or mopping with water for injections**

This relates partly to the techniques used and partly to the quality of the cleaning materials. With the techniques for the cleaning and disinfection of cleanroom floors, either the "two-bucket" or the "three-bucket" technique is recommended. Both of these techniques involve using a bucket of disinfectant and a bucket of water.

In the "two-bucket" technique there is a "wringer" (for the mop) over the bucket of water. In the "three-bucket" technique there is a third bucket, empty except for having a wringer mounted over it. To illustrate this, the three-bucket technique involves:

1. Dipping the mop in disinfectant
2. Mopping the floor
3. Dipping the mop in the bucket of water
4. Rinsing the excess water off the mop head into the third (empty) bucket
5. Dipping the mop in disinfectant
6. Repeating

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4.4 Deeper cleaning

Sometimes a major decontamination exercise is required for a cleanroom (such as following an HEPA filter failure or following maintenance work). There are different techniques used to decontaminate cleanrooms. Traditionally this has involved either manual cleaning or disinfection by operators using mops and buckets, or, in more recent years, decontamination with fumigation units and chemicals such as hydrogen peroxide or peracetic acid (noting appropriate health and safety regulations). For this method, gas generators can be placed in cleanrooms at known locations. This takes some time and effort in mapping the cleanroom (which additionally requires some form of validation using chemical and biological indicators).

To aid such techniques commercial bucket and trolley systems are available. For aseptic filling areas it is important that these are manufactured from a material that can be autoclaved, such as stainless steel. In terms of the materials used for cleaning, appropriate mops and wipes should be used. Wipes should be sterile and be low particle shedding. When using wipes for surfaces, the detergent or disinfectant should ideally be sprayed onto the surface rather than the wipe and left for the required contact time before wiping is undertaken.

The design of the mop is a further consideration and must be easy for the cleanroom operator to use and allow for a uniform application of the disinfectant solution.

4.5 Cleaning procedure of production rooms

Some of the procedures for cleaning of equipment and production rooms will be in the following sequence:

I Cleaning Frequency

1. Cleaning activities are performed before initiating

a manufacturing cycle and following the completion of a manufactured batch of product.

2. Where production areas are not in use for extended periods, cleaning steps must be performed at least every sixty (60) days.
3. Rooms used in equipment cleaning operations are to be cleaned biweekly (every two weeks) or more often as activity/use dictates.
4. Any storage rooms or office areas that are located within a production area must also be cleaned weekly. Floors must be vacuumed or swept and mopped. Surfaces, table tops, cabinet tops etc. must be wiped with Simple Green (3 oz. /gallon water).
5. Cleaning activities may be performed more frequently as determined necessary by the good judgment of operations and production personnel.
6. If maintenance activities compromise cleaning the area, supervisor will schedule additional cleaning as needed.
7. Forms (specific documents/records) will identify the date of area cleaning activities

Note: Room cleaning logs will be maintained for individual rooms, and in some cases by area or production line, where single activities are performed in multiple rooms or area.

II Pre-cleaning Activities, Where Applicable

- Remove all materials (drums, components, etc.) remaining in the area after the completion of a manufacturing cycle for proper disposition.
- Remove portable equipment (portable pumps, frames, etc.) from the area after the completion of a manufacturing cycle for proper disposition/cleaning.

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- Remove rubber mats from the area.
- Sweep or vacuum the floors as appropriate for the area.
- Where necessary, cover electrical outlets and turn power switches on production equipment off.

Learning activities

1. Research further on cleaning procedures for Cleanroom
2. Demonstrate cleaning cleanrooms in real life pharmaceutical setting (Simulated Exercise)

Test yourself

What are the major factors that need to be considered while sterilising cleanrooms in pharmaceutical firms?

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Unit 5

Documentation of the Cleaning Process

Learning outcomes

When you complete this Unit, you will:

1. Understand the importance of maintaining of records of cleaning and validation.
2. Apply cGMP procedures of cleaning process.

5.1 Introduction

The regulatory authorities expect firms to have written procedures (SOP's) detailing the cleaning processes used for various pieces of equipment. If firms have one cleaning process for cleaning between different batches of the same product and use a different process for cleaning between product changes, they expect the written procedures to address these different scenarios. Similarly, if firms have one process for removing water soluble residues and another process for non-water soluble residues, the written procedure should address both scenarios and make it clear when a given procedure is to be followed.

5.2 Written cleaning process

Depending upon the complexity of the system and cleaning process and the ability and training of operators, the amount of documentation necessary for executing various cleaning steps or procedures will vary.

When more complex cleaning procedures are required, it is important to document the critical cleaning steps (for example certain bulk drug

synthesis processes). In this regard, specific documentation on the equipment itself which includes information about who cleaned it and when is valuable. However, for relatively simple cleaning operations, the mere documentation that the overall cleaning process was performed might be sufficient.

Other factors such as history of cleaning, residue levels found after cleaning, and variability of test results may also dictate the amount of documentation required. For example, when variable residue levels are detected following cleaning,

The Arab good manufacturing practice guidelines (1995)

5.3.6.10 b.A procedure for cleaning production equipment should specify cleaning of the equipment prior to change of a batch as well as change of a product. The procedure should include the method of cleaning, the tools and cleaning materials to be used. The cleaning operation should be documented.

5.3.6.10 c.A procedure for cleaning of manufacturing areas should be available. This procedure should include the specific areas to be cleaned, the tools and cleaning materials to be used, the time and sheue of cleaning. **The cleaning operation should be documented.**

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particularly for a process that is believed to be acceptable, one must establish the effectiveness of the process and operator performance. Appropriate evaluations must be made and when operator performance is deemed a problem, more extensive documentation (guidance) and training may be required.

5.3 Cleaning records

cGMP guidelines demand that a written record of major equipment cleaning, maintenance (except routine maintenance such as lubrication or minor adjustments) should be include in individual equipment logs that show the date, time, product, and lot number of each batch processed.

In case where dedicated equipment is employed, the records of cleaning, maintenance and use shall be part of the batch record.

The persons performing and double checking the cleaning and maintenance should date and sign the log indicating the work performed. Entries in the log shall be in chronological order.

Cleaning activities are to be recorded in the appropriate cleaning log and/or on the applicable cleaning activity form. Each entry will be initialed and dated by the individual performing the cleaning activity.


Fig 5.1 Equipment Cleaning Details (e-version)

A second individual must verify completion of the cleaning activity by initial and date.

All manufacturing rooms and equipment (except stainless steel utensils and scoops, spatulas, etc.) have a cleaning record which documents products processed in the room/equipment, batch number, cleaned by and checked by, and date.

JORDAN PHARMACEUTICALS Pvt Ltd. EQUIPMENT CLEANING RECORD

Equipment No	P0000			Equipment Type	XYS	
Cleaning SOP Ref	000000			Issue No	P0000	
Date	Cleaned by	Previous product	Batch No	Type of Clean	Initial	Date

Table 5.2. Equipment Cleaning Record Format

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5.4 Labelling of cleaning agents

All cleaning agents and disinfectants shall be labelled with 'Not for Pharmaceutical Use' labels. See below:

NOT FOR PHARMACEUTICAL USE

FOR CLEANING ONLY NOT A PHARMACEUTICAL INGREDIENT

The Production Operator applying the label should complete all fields on the label, or insert N/A, if not applicable. No blank fields shall be left on the label.

Whenever there is equipment failure and the equipment is removed from a processing room intentionally (unless being transferred from one processing room to another for continuation of processing), the operator must inspect the equipment and remove any granules/tablets/capsules and place them into an appropriate form of holding container/bucket/drum.

If out of specification or contaminated product is found then it should be placed into a controlled waste bucket and the existing clean label should be removed

from the equipment and adhered to the batch manufacturing record.

The details of the previous product should be entered into the relevant clean record log book and the operator is to ensure that one quarantine hold label remains on the equipment prior to removal of equipment from the room.

Learning activities

Demonstration on filling up of cleaning and sanitization records.

Exposure to the reporting procedure related to cleaning.

Test yourself

Explain the process in cleaning and sanitization.

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