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VOCATIONAL TRAINING CORPORATION (VTC)

**TECHNICAL OPERATIONS IN THE PHARMACEUTICAL
PRODUCTION UNITS
(LEVEL 3)**

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Acknowledgment

This publication has been made possible with the invaluable contribution by many stakeholders including competitive inputs and leadership by the pharmaceutical industry in Jordan, instructors, and consultants.

The Ministry of Labour and ETVET Council through Vocational Training Corporation (VTC) and the European Union with technical assistance from GOPA Worldwide Consultants have partnered with the following industry experts, the Jordanian Association for Pharmaceutical Manufacturers (JAPM) to develop the curriculum framework and the resulting learning materials.

The VTC would like to thank particularly the following Consultants and Members of the Curriculum Development Panel for their works in research, developing the content, review and finalising these materials.

Dr. Hanan Sboul	JAPM
Mr. Rami Dari	Ram Pharma
Eng. Sana' Abbas	JPM
Dr. Ayyad Rumman	Ayyad Rumman Pharma Consulting Centre
Mrs. Suzanne Grigoleit	VTC
Dr. Judy Saleh	PCOE - VTC
Eng. Hussain Zuhairi	PCOE - VTC
George Kolath	Lead ISTE, GOPA

No part of this publication may be reproduced, stored in a retrieval system, or transmitted in any form or by any means, electronic, mechanical photocopy and recording or otherwise, without the written permission of Vocational Training Corporation.

The publishers have made every effort, where required, to acknowledge copyright holders, but if any have been overlooked, the necessary arrangements will be made at the first opportunity. We welcome any information that will enable us to correct, in subsequent edition, any inaccurate or omitted credit.

No liability is assumed with respect to the information contained in.

Published by

Vocational Training Corporation
Amman - Jordan

With financial and technical assistance by the European Union and GOPA Worldwide Consultants respectively.

© Vocational Training Corporation
2014

All rights reserved.

This Project is Funded by the European Union

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Overview of the Module

The purpose of this Module is to explain processes and procedures in encapsulation.

This Module, when completed, will provide you with the necessary competencies to apply theoretical and practical concepts in the formulation process of encapsulation by applying appropriate principles and procedures. They include the correct and safe handling, used and cleaning of different types of encapsulation machines and correct monitoring the machines and required procedures during the production process.


Learning Outcomes

When you complete the Module, you will have the following specific competencies:

1. Perform encapsulation through the application of encapsulation
2. Apply correct and safe handling, use and cleaning of different types of encapsulation machines
3. Employ correct documentation procedures and apply SOPs correctly.

Credit Value

The credit value of this Module is 5 credits which equals to 105 hours of nominal learning time. This includes theory lesson (face-to-face), supervised and self-directed practice sessions, and continuous and summative assessment (learning controls).

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Assessment

At the end of every Unit there are learning and self-assessment activities. The learning and assessment activities provide you with the opportunity to apply knowledge and skills you have acquired. You should successfully complete all the learning and assessment activities as evidence that you actually master the skills and related knowledge in this module.

Your Instructor too will carry out a formal assessment as evidence of your certification. As before, the mode of assessment will be both continuous (formative assessment) and summative (external assessment).

In this Module, your continuous assessment, both theoretical and practical consisting of 30% of all units, will be considered for the final grading.

The final assessment of the entire Module will be an examination of both theory (30%) and demonstration or practicals (40%), which will be carried forward for the purpose of certification.

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Learning Outcomes

At the end of this Unit, you will:

1. Interpret the physical and chemical characteristics of capsules and their significance in encapsulation.
2. Read and apply the SOPs on encapsulation.

1.1 Introduction

The encapsulation process is a procedure where a machine, which is fed with a blend of ingredients and empty hard gelatin capsules, opens the capsules, fills them to a set weight with the ingredient powder and closes the capsules again. Different machines work in slightly different ways and with differing levels of automation.

Capsules are popular solid dosage form since they have a certain elegance (colourful and appealing), ease of use, and portability.

In this Unit, we analyse the principles of encapsulation including their physical and chemical properties.

1.2 Definition

The word 'capsule' is derived from the Latin *capsul*, which means a *small box*.

In the pharmaceutical industry, capsule is used to describe an edible package made from gelatin or other suitable material which is filled with medicines

Mothes and **Dublanc**, two Frenchmen, are credited with the invention of the gelatin capsule, patented in **1834**.

In **1848**, **James Murdock** developed the two-piece telescoping capsule, as seen today.

to produce a unit dosage.

In other words, *Capsules are solid dosage form in which the drug substance is enclosed within a hard or soft soluble shell generally formed from gelatine.*

Two types of capsules are available in the industry today: **hard** and **soft gelatin capsules**.

1.3 Advantages of Capsules

The capsules have some advantages over other dosage forms and let us state them:

1. Elegance, ease of use and portability.
2. They provide smooth, slippery, tasteless shell to swallow.
3. They are convenient method by which liquids may be orally administered to patients as a unit dosage form.
4. Drugs having unpleasant taste and odour can be administered through capsules.

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5. They have better bioavailability than tablets.
6. Possible to design sustained or prolonged release formulation in capsule form.
7. Economical (low cost).
8. They are difficult to be counterfeited.
9. The stability of therapeutic agents may be improved in a capsule formation.
10. They are available in various shapes and sizes, and allow flexibility in formulations.

1.4 Disadvantages of Capsules

1. Extremely soluble drugs (e.g. potassium iodide or potassium bromide) cannot be administered in capsules as its sudden release can cause irritation/burning sensation.
2. Problems associated with homogeneity of fill weight and content may be associated with capsule formulations.
3. The requirement for specialised manufacturing equipment.
4. Potential stability problems associated with capsules containing liquids.

1.5 Types of Capsules

Capsules are classified as hard gelatin capsules (HGC) and soft gelatin capsules (SGC). Capsules are also categorised as 'two-piece' or 'single piece' capsules.

The hard capsule consists of two pieces in the form of cylinders closed at one end; the shorter piece, called the 'cap' fits over the open end of the longer piece, called the 'body'.

Hard-shelled capsules are typically made using gelatin and contain dry, powdered ingredients or miniature pellets (e.g. processes of extrusion or spherulisation).

Soft-shelled capsules are primarily used for oils and for active ingredients that are dissolved or suspended in oil. The soft capsules are also known as 'soft gels'.

1.6 Differences between Hard and Soft Gelatin Capsules

Table 1.1 Hard Gelatin versus Soft Gelatin Capsules

Hard Gelatin Capsules		Soft Gelatin Capsules
a	Shell does not contain plasticiser	Shell contains plasticiser
b	Available mostly in tubular or oblong forms	Available in any shapes
c	They are not hermetically sealed and may leak at the joints	Hermitically sealed, no leakage possible
d	Ease of production on small scale	Complex production methods/machines and so large scale is economical
e	Mostly for oral, rectal or vaginal insertion	They are multi-functional- oral, topical, ear nose, rectal and vaginal insertion
f	Meant for enclosing solid medications	Flexible for both liquid and semisolid medications
g	Construction is two-step process where shell is made on one type of machine and filling is on another type of machine	Construction, filling and sealing are done on single machine
h	Accuracy of the fill content is less as compared to the soft gelatin capsules, as powder form is complex.	Since fill content is liquid or semi-solid, it is accurate and precise filling is possible

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
1.7 Physical and Chemical Properties of Capsules

Capsule Size	Volume (ml)	Mg of Lactose	Mg of Aspirin
000	1.37	1340	1000
00	0.95	929	600
0	0.68	665	500
1	0.5	489	300
2	0.37	362	250
3	0.3	293	200
4	0.2	195	125
5	0.13	127	60

Similar raw materials have been used in the manufacture of both types of capsules. Traditionally both contain **gelatin, water, colorants** and optional materials such as **process aids and preservatives**. The soft capsule contains various plasticizers, such as glycerine and sorbitol.

Mostly today, hard capsules are produced from

hypromellose (in order to produce a shell with low moisture content) and soft capsules are made from *modified starch*.

1.7.1 Gelatin

Gelatin is a major component used for capsules. Gelatin possesses five basic properties that make it suitable for the manufacture of capsules:

1. It is non-toxic, widely used in foodstuffs, and acceptable for use throughout the world.
2. It is readily soluble in biological fluids at body temperature.
3. It is a good film-forming material, producing a strong flexible film. The wall thickness of a hard gelatin capsule is about 100 μm .

4. Solutions of high concentration, 40% w/v, are mobile at 50 C.
5. A solution in water or in a water-plasticizer blend undergoes reversible change from a sol to a gel at temperatures only a few degrees above ambient. This is in contrast to other films formed on dosage forms where either volatile solvents or large quantities of heat are required to cause this change of state (e.g, tablet film coating).

Gelatin is a substance of natural origin. It is prepared by hydrolysis of collagen, which is the main protein constituent of connective tissues.

Animal skins and bones are the primary raw materials used for the manufacture of gelatin.

There are two types of gelatin that are popular in the industry:

- a) produced by acid hydrolysis
- b) produced by basic hydrolysis.

The acid process takes about 7-10 days. The basic process takes about 10 times as long as acid process since basic process is used mainly for bovine bones.

After hydrolysis, the gelatin is extracted from the treated material using hot water. The first extracts contain the gelatin with the highest physical properties and as the temperature is raised, the quality falls. The resulting weak solution of gelatin is concentrated in a series of evaporators and then chilled to form a gel. This gel is then extruded to form strands, which are then dried in a fluidised-bed system. The dried material is graded and then blended to meet the various specification requirements.

The properties of gelatin that are most important to the capsule manufacture are the Bloom Strength and Viscosity. The Bloom Strength is a measure of

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gel rigidity. It is determined by preparing a standard gel and maturing it at 10 C. It is defined as the load in grams required to push a standard plunger 4 mm into the gel.

1.7.2 Colourants

The colourants that can be used are of two types: water soluble dyes or insoluble pigments.

To make a range of colours, dyes and pigments are mixed together as solutions or suspensions. The dyes are mostly synthetic in origin and can be subdivided into the azo dyes and non-azo dyes.

Two types of pigments are used: a) iron oxides, black red and yellow; and, b) titanium oxide, which is white and used to make the capsule opaque.

1.7.3 Process Aids

Sodium lauryl is used as a process aid in hard capsules as a wetting agent to ensure that the lubricated metal moulds are uniformly covered when dipped into the gelatin solution.

Preservatives are added to the hard capsules in order to prevent microbiological contamination during manufacture.

1.8 Soft Gelatin Capsules

The Softgel is a hermetically sealed, one-piece capsule shell with a liquid or semisolid fill without a bubble of air or gas. The fill materials can include a wide variety of vehicles and can be either a solution or a suspension. Though the Softgels may be either clear or opaque, it is standard practice to use a clear shell (clear coloured or natural amber) only when the fill is also a clear solution.

Advantages of hard gelatine capsules		over soft gelatine capsules	
Contain 4-5 times less gelatine than soft gelatine capsules			Require 4-5 times more gelatine than the hard gelatine capsules
Require no other additives. Consists of water and gelatine only			Require addition of glycerin for softening purposes
Allow step-by-step filling of 2 different formulations (i.e. 2-stage-release)			Have to be sealed immediately after filling one substance (filling and sealing are one and the same process)
Heat resistant : allow filling of thermo-stable substances up to 75°C			Filling temperature limited to about 35°C : filling of solid substances with higher melting points impossible
Are stable in hot climates			Tend to stick together and become gluey
Will disintegrate faster due to the capsule wall being five times thinner than the walls of soft gelatine capsules			Will disintegrate slower due to the thickness of its gelatine/glycerin wall
Less product migration into the shell, less diffusion of odours			Glycerin acts as a plasticiser by disrupting the gelatine structure - consequently, higher diffusion into and through the walls
Constant external dimensions (easier blistering/packaging)			Dimensions vary according to filling weight and vary throughout a batch

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Advantages of hard gelatine capsules		over soft gelatine capsules
Contain 4-5 times less gelatine than soft gelatine capsules		Require 4-5 times more gelatine than the hard gelatine capsules
Require no other additives. Consists of water and gelatine only		Require addition of glycerin for softening purposes
Allow step-by-step filling of 2 different formulations (i.e. 2-stage-release)		Have to be sealed immediately after filling one substance (filling and sealing are one and the same process)
Heat resistant : allow filling of thermo-stable substances up to 75°C		Filling temperature limited to about 35°C : filling of solid substances with higher melting points impossible
Are stable in hot climates		Tend to stick together and become gluey
Will disintegrate faster due to the capsule wall being five times thinner than the walls of soft gelatine capsules		Will disintegrate slower due to the thickness of its gelatine/glycerin wall
Less product migration into the shell, less diffusion of odours		Glycerin acts as a plasticiser by disrupting the gelatine structure - consequently, higher diffusion into and through the walls
Constant external dimensions (easier blistering/packaging)		Dimensions vary according to filling weight and vary throughout a batch

The standard Softgel shapes for oral pharmaceutical products are oval, oblong, and round. The size of the Softgel is represented by a numerical number, which represents its nominal capacity in minims (1 cc j 16.23 minims). For example, an 11 oblong Softgel can be filled with 8.5–11.0 minims of fill formulation.

Softgels can be easily manufactured in any shape with a plane of symmetry and any size (to contain up to 25 ml) via appropriate die design. The shapes of soft gelatin capsule are round, oval, oblong, tube. A recent survey has shown that smaller sized Softgels

are preferred within each shape category, with oval being the most popular shape.

1.9 Application of Soft Gelatin Capsules

The pharmaceutical applications of soft gelatin capsules are:

- Oral dosage form
- Suppository
- Specialty Packaging for single use application such for ophthalmic, and ointment products


Fig 1.1. Softgel components

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

Fig 1.2 Examples of softgels

1.9.1 Advantages of soft gelatin capsules

Some of the distinctive advantages of soft gelatin capsules are:

- They permit liquid medication to become easily portable.
- Accurate, precise and uniform dose.
- Hermetically sealed (inherently)
- Better bioavailability of the drug and fast onset of action
- Biopharmaceutical characteristics can be adjusted by selecting liquid, emulsion, suspension or powder fill.
- Reduced dustiness; lack of compression stage in manufacture
- Reduced gastric irritancy compared to tablets and hard shell capsules
- Safe: reduced operator and environmental contamination with highly potent or cytotoxic compounds
- Smooth slippery and easy to swallow

- Specialty packaging become available

1.9.2 Disadvantages of Softgel Capsules

- Generally, product is contracted out to a limited number of specialty houses, e.g. Scherer, Banner.
- Generally more costly to produce than tablets or hard shell capsules
- More intimate contact between the shell and contents than with dry-filled hard shell capsules - stability a concern.
- Not adaptable to incorporation of more than one kind of fill into the same capsule (compare with hard shell capsules)

ROUNDS						
Height: Diameter: CC	2 1.2-2.0 0.072-0.174	3 2.2-3.0 0.102-0.188	4 2.0-3.0 0.172-0.244	5 2.2-3.0 0.107-0.208	6 4.0-6.0 0.281-0.330	7 5.0-7.0 0.508-0.874
OBLONGS						
Height: Diameter: CC	2 2.5-3.0 0.142-0.165	3 2.0-4.0 0.105-0.244	4 4.0-6.0 0.244-0.309	5 5.0-6.0 0.305-0.370	6 8.0-9.0 0.405-0.493	11 8.0-11.0 0.824-0.870
OVALS						
Height: Diameter: CC	2 1.5-1.8 0.097-0.111	3 2.4-3.0 0.162-0.185	4 2.4-3.0 0.174-0.249	5 4.2-5.0 0.328-0.309	6 5.2-6.0 0.523-0.570	7 7.44 0.302-0.342
TUBES						
Height: Diameter: CC	2 2.4-3.0 0.142-0.309	3 5.0-6.0 0.308-0.370	4 5.0-6.0 0.375-0.493	5 11-17 0.722-1.047	11 17-20 1.463-1.571	20 22-25 1.463-1.571
SUPPOSITORIES						
Height: Diameter: CC	2 5.0-6.0 0.508-0.570	10 7.0-10.0 0.450-0.620	20 10.0-15.0 1.000-2.400	60 22.0-25.0 3.817-4.500		

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1.10 A Comparison between Tablet and Capsule Production Stages

In comparison with conventional tablet production, capsule production does not require expensive and time consuming operations like repeated mixing and sieving, or granulation and compression. This makes capsules as an attractive option for the pharmaceutical industry.

Table 1.2 Production process operations in the manufacture of tablets and capsules.

Tablets	Capsules
<ol style="list-style-type: none">1. Weighing2. Preparing ingredients3. Mixing4. Granulating5. Drying6. Sieving7. Addition of lubricants/mixing/sieving, as necessary8. Compression9. Packing	<ol style="list-style-type: none">1. Weighing2. Preparing ingredients3. Mixing4. Filling into capsules5. Packing

Learning Activities

Carry out further internet research on the principles and emerging technologies in encapsulation.

Test Yourself

1. What are the principles and functions of encapsulation?
2. Define the physical and chemical properties of capsules.

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Learning Outcomes

At the end of this Unit, you will:

1. Illustrate the procedures for encapsulation by interpreting SOPs.
2. Feed the capsulation machine with the mixed product.
3. Assemble and demonstrate the correct encapsulation procedures.

2.1 Introduction

You have learnt already the principles of encapsulation and are able to apply the SOPs for encapsulation. In this Unit, you will carry out the process of encapsulation on different encapsulation machines.

2.2 Encapsulation Machines

Commonly referred to as **capsule filler**, the encapsulation machine has the ability to fill many different products. Powders, granulations, liquids, tablets and capsules can be filled into a two piece capsule.

Encapsulation machinery technology varies a great deal from one manufacturing to the next. Not all machines can fill a wide variety of products; most are designed to handle free flowing powders much like powders that are prepared for a tablet press.

The capsule filler must first position all of the incoming capsules into an upright position

(rectification), separate the cap from the body (top from bottom), attain the proper fill volume (capsule weight), and then the product filled body is rejoined with the cap and ejected from the machine. Some capsule filling machines have the ability to compress or tamp the powder for proper filling volume and weigh control.

Encapsulators can be defined as:

- 1) Hand operated
- 2) Semi-Automatic
- 3) Fully Automatic

The Hand Operated capsule filler requires the operator to organize the capsules in the correct position, separate the cap from the body, and fill and close the caps (basically the hand filler is a holder for the capsule body).

There are exceptions and some hand operated fillers assist the operator with separation and closing functions. The Semi-automatic machine requires the operator to move rings (capsule holder rings) from the rectifier to the filling and closing stations allowing for production up to 25,000 capsule per hour.

Automatic machines with speeds up to 90,000 per hour can be divided into two categories: Continuous and Intermittent operations. The intermittent motion machine is divided into segments. Each segment indexes from each machine function; rectify, fill, tamp, close and eject.

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The automatic machine is a continuous operation somewhat comparable to a rotary tablet press in that the rotation is continuous and does not start and stop.

2.3 Capsule Filling

The dosing systems in industrial scale filling machines may be divided into two groups:

I Dependent (Auger) systems: which use the capsule body directly to measure the powder. Uniformity of the fill weight can be achieved if the capsule is filled completely.

The dependent dosing systems are semiautomatic in operation, requiring an operator to transfer the capsule holder from one operation to the next. The output from such machines varies between 15000 and 25000 per hour and is dependent on the skill of the operator.

II Independent systems: where the powder is measured independently of the body in a special measuring device. Weight uniformity is not dependant on filling the body completely and the

capsule can be part filled.

The independent dosing systems are fully automatic systems, using a dosing mechanism that forms a 'plug' of powder that is transferred then to the capsule shell.

Plugs are soft compacts formed at low compression forces (10-100N), they are soft because they are not the final dosage form as they will be contained inside the capsule shell.

There are two types of plug forming machines:

- Dosator system:
- Tamping finger and dosing disc system.

Dosator systems:

These consist of dosing tube inside which there is a movable spring-loaded piston, thus forming a variable- volume chamber in the bottom of the cylinder.

The dosator system is the most commonly used capsule filling system and could either be intermittent or continuous in motion. The output could reach up to 150 000/ hour .


Fig 2.1 Diagram of a dosator or dosing tube system
(A) compression force platen;(B) piston; (C) dosing tube; (D) powder hopper; (E) plug ejection platen;(F) capsule body in bush; and (G) powder plug.

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
2.4 Tamping Finger and Dosing Disc

- In the system a dosing disc forms the bottom of a revolving powder hopper.
- The dosing disc has a series of accurately drilled holes in which powder plugs are formed by several sets of tamping fingers.
- Sets of tamping fingers (stainless steel rods) are lowered into them through the powder bed to compress the plug.
- At each position the fingers push materials into the holes, building up a plug before they index on to the next position.
- Thus, the plug is formed in a series of tamps and not in a single action as on the dosator machines
- At the last position the finger pushes the plug through the disc into a capsule body.
- The powder fill weight can be varied by the amount of insertion of the fingers into the disc, by changing the thickness of the dosing disc, and by adjusting the amount of powder in the hopper. Thus selection of the correct thickness of disc is important, because if the fill weight cannot be achieved, the machine has to be dismantled to change it
- The machines that use this system are all intermittent in motion.

2.5 Pellet Filling

Capsules for modified release purposes are usually formulated with granules or coated pellets fillings. They are filled on an industrial scale using machines adapted from powder use (pellets are not to be compressed during the process). All dosing systems are based on chamber with the volume that can be easily changed (Dosator systems). Since pellets are not supposed to be compressed, they are held inside the measuring devices by mechanical means


Fig 2.2 Schematic diagram of a dosing disc and tamping finger system:

(A) over-load relief spring; (B) tamp-depth adjuster; (C) tamping finger; (D) powder hopper; (E) powder bed; (F) dosing disc; (G) suction; (H) support plate; (I) ejection adjuster; (J) guide block; (K) transfer block; (L) capsule body in bushing; and (M) powder plug.

(e.g. inverting the dosator or applying suction to the dosing tube).

In calculating the weight of pellets that can be filled into a capsule, it is necessary to make an allowance for their size. Unlike powders (which have a much smaller size), they cannot fill all available space within the capsule because of packing restrictions.

The degree of this effect will be greater the smaller the capsule size and the larger the particle diameter.

2.6 Procedures for Preparation of Raw Materials for Hard Gelatin Capsules

The first step in the process is the preparation of the raw materials.

- A concentrated solution of gelatin, 35-40%, is prepared using demineralised hot water (at 60-70 C), in jacketed pressure vessels.

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- This solution is stirred until the gelatin has dissolved and then a vacuum is applied to remove any entrapped air bubbles.
- Aliquots of this solution are then dispensed into suitable container and required amount of dye solutions and pigments suspensions added.
- The viscosity is measured and adjusted to a target value by the addition of hot water (remember, higher the viscosity, thicker the shell wall to be produced).
- The prepared mixes are then transferred to a heated holding hopper on the manufacturing machine.
- Hypromellose solutions are prepared in a similar way but with the additional gelling agent such as carrageenan, and a co-gelling agent, such as potassium chloride, to convert them into a gelling system.

2.7 (A Model) SOPs for the Encapsulation Process

1. The Manufacturing room must be cleaned and signed off by a QA officer, Production Supervisor (PS) or Team leader (TL).
2. The cleaning of the encapsulation machine must be approved prior to and after the installation of the change parts and ancillary machine parts.
3. The status of the following ancillary equipment for the encapsulation machine must be approved by a QA Officer, PS or TL before assembly on the machine can take place:
 - Powder filling station.
 - Powder/granulate supply hopper.
 - Capsule supply hopper
 - Faulty capsule station.
 - Closing station
 - Capsule discharge station.

- Extraction hoses
- Compressed air line(s)

The equipment must be left in the clean encapsulation room on the appropriate tool trolley provided ready for inspection and approval.

All relevant labels must be approved and signed for prior to use of the equipment.

4. The product-specific 'Capsule Size Change Parts' must be checked and approved by a QC Officer, PS or TL before they can be brought into the room. These parts comprise the following:
 - Tamping Pins
 - Capsule segments (upper and lower)
 - Dosing Disc
 - Capsule Magazines
5. Place a cleaned, approved discharge cone and lay-flat tubing in the space provided in the corresponding discharge station over the encapsulation room.

Ensure that the tubing is secured to the base of the discharge cone and the powder/granulate supply hopper.

Place an 'IN USE' status label for the discharge station and a separate 'IN USE' status label for the discharge cone on the discharge control panel.

The 'In use' label must indicate:

- Product name
- Product strength
- Product batch number
- Part / room number

6. Bring the following into the Manufacturing room:
 - The appropriate BMR.
 - Clean buckets, labels, clear polythene and black bucket liners.
 - Labels (Bucket 'IN-PROCESS' and 'IN USE')

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- Sample containers
 - Polypropylene scoop.
 - Floor balance
 - Table balance
 - Capsule transfer system
 - Capsule de-duster/ metal detector and sorter (if applicable)
 - Capsule polisher (if applicable)
7. Insert the required liners into the buckets and attach the bucket labels. Unless otherwise advised by the PS/TL or QC officer, each bucket will require one black liner, one clear liner.
 8. The QC Officer, PS or TL must check the following:
 - All equipment and materials are clean.
 - All information is correctly recorded.
 - All details of materials within the room are entered into the BMR.
 - All buckets are checked as being clean prior to placing liners inside buckets.
 9. A note must be made in the BMR regarding the total number of approved buckets at the start of encapsulation. If for any reason during encapsulation additional buckets are required, these buckets are to be independently checked prior to use and a comment inserted into the relevant section of the BMR.
 10. Locate the Intermediate Bulk Container (IBC) containing the granulate/powder that is to be encapsulated and move it to the correct drop station. Do not load the IBC onto the discharge station at this time. It is necessary to have the identity of the granulate/powder and the IBC confirmed by a QC officer, PS or TL before the IBC can be loaded onto the discharge station.
 11. Load the IBC onto the discharge station using the lifting equipment. The IBC must have also an 'IN USE' status label attached indicating the product name, strength, batch number and IBC number.
 12. Place the "IN-USE" labels on the white buckets.
 13. Place 'IN USE' status labels on the following:
 - Door of the Manufacturing room
 - Encapsulation machine
 - Polisher/de-duster, metal detector and sorter
 - Table balance
 - Floor balance
 - Vacuum cleaner
 - Drop station
 - Drop cone
 - Capsule conveyor
 - Capsule polisher
 14. The batch must be signed into the room by an QA officer, PS or TL before encapsulation can begin.
 15. Connect the compressed air at the service assembly.
 16. Charge the capsule supply hopper with empty capsule shells.
 17. Charge the Powder/granulate supply hopper with material from the IBC.

2.8 Risks and Hazards Associated with the Encapsulation Processes

In Module 102, we discussed extensively on various occupational health and safety measures to ensure that you are well protected to limit your exposure to various safety risks and hazards during various manufacturing operations. In the encapsulation process, like the tableting process, you are faced with many similar safety risks and you are advised to take appropriate caution to protect yourselves.

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Remember at all stages, that appropriate PPE should be worn, safety precautions adhered to, and that you are constantly monitoring your health through regular health check-ups too.

Some of the following risks are associated with the encapsulation process:

- Fast-moving parts, high noise levels and dust emissions during their operation of the machines/equipment.
- Special concerns when workers are exposed to toxic solvent vapours and potent drugs as airborne dusts during the unit operations of capsule filling.
- Capsule products generate solid and liquid wastes during cleaning and sterilization, and

from leaks and spills and rejected products.

- Waste waters may contain inorganic salts, sugars, and traces of drug substances or powders which may cause both environmental pollution and safety risks to the workers.

2.9 Significance of Monitoring the Encapsulation Machine

The monitoring and continuous quality checks on the encapsulation machine will ensure higher quality capsules with no defects. This is also carried in order to ensure that the machine is operating smoothly without any faults.

The following table highlights some of the

2.9.1 Manufacturing problems

FAULT	CAUSE	MACHINE REMEDY
Machine stops	'Sticky' product	Optimise height configuration of the dosing disk
		Tamping ring facing dosing disk at an angle
		Provide recesses at the dosing disk
		Increase air gap (S~0.1mm)
De-mixing of powder	Powder has bad flow properties	Optimise tamping pressures.
		Optimise powder feed to the tamping stations for example by additional powder levellers (to be provided by Engineering)

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Powder is ejected from the dosing bore during tamping	Powder with air enclosures	Adjust the powder level
		Adjust the number of cycles
		Adjustment of tamping pressures: Maximum tamping pressure at tamping station Lower tamping pressure at the following tamping stations ->a>b>c>d>e
		Review the possible use of softer tamping pin springs
Powder collapses upon ejection from dosing bore	Powder with a dust percentage	Air gap (S, 0.02mm)
		Optimise tamping pressures
		Reduce gap (S~0.02-0.03mm)
		Review the possible use of 'Harder' tamping springs

faults, causes and how to remedy them on the encapsulation machine.

You should also be aware that various quality assurance tests are carried out in order to ensure quality finished products. This is also true in the case of encapsulation processes where the Production Team should carry out 'start up analysis tests' and 'end analysis and finish-up tests'. These checks can be categorised as IPCs too.

2.9.2 Start-Up Analysis Tests

The following aspects are analysed at the start of the encapsulation process:

a. Capsule weight

1.1 Weigh 20 capsules individually using the IPC checkweigher located in the In-Process Control Room (IPC).

1.2 Review average filled capsule weight together with high and low weights.

Calculate the filled capsule weight bias on the IPC by summing the high and low weight deviations from

the average.

1.3 Write the result on the IPC printout. Ensure that the batch number, product code, 'T' limits (T1, T2 and T3) and target weights are correct on the IPC printout as per the BMR.

1.4 Carefully review the 'flag' weights from the printout against the BMR specifications. The individual filled capsule weights must comply with the limits as specified in the BMR.

1.5 Attach the printout to the BMR 'Start-Up Analysis' page.

b. Capsule length

2.1 Measure the length of 10 capsules with a digital callipers or suitable measuring equipment.

2.2 Record each individual result and the mean value on the 'Start up Analysis' page.

c. Capsule Disintegration

3.1 Perform a disintegration test on 6 capsules.

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3.2 Record the results on the 'Start-Up Analysis' page.

d. Capsule condition check

Check the overall condition of an appropriate sample size of capsules for any of the following defects (the sample size is defined on the BMR):

- Wrong capsules
- Empty capsules
- Wrong capsule contents
- Foreign matter present
- Split/Broken capsules
- Leaking capsules
- Underweight capsules
- Marked/Dented capsules
- Capsule Closure Failure
- Poorly polished capsules

Record the results in the 'Tablet Condition Check' section of the 'Start up Analysis' page.

Ensure that all the test results that are entered into the BMR at the product 'Start-Up Analysis' and throughout the encapsulation section are rounded correctly to the required number of decimal points, in comparison to the BMR specification.

Note: If an Out-of-Specification (OOS) result is encountered during the 'Start up Analysis', the machine may require adjustment to correct the OOS result.

Ensure that the correct numbers of capsule samples are taken for both the 'Start-Up Analysis' and throughout the encapsulation run.

A QA officer, PS or TL must approve and sign off the 'Start Up Analysis' sheet before production can begin, to confirm that the capsules are being manufactured to the BMR specifications.

All capsules produced up to this point must be

discarded into a suitable waste container, which must be labelled accordingly. All samples weighed on the table balance must be disposed of in the 'Capsule Reject Bin' located on the wall beside the table balance. Unless otherwise requested, the 'Start up Analysis' sample must be disposed of in the Parma waste bin in the IPC.

After completion and approval of the 'Start Up Analysis', set up the ancillary equipment to the front of the machine so that the capsules can undergo de-dusting, metal checking, sorting and polishing (if required) are de-dusted on line. Once the metal checker is installed the operator must carry out a metal token challenge. This involves passing a 0.3mm FE, 0.5mm Non-Fe and a 0.5mm SS token through the metal checker. This must be recorded in the Manufacturing BMR. (If the challenge fails the operator must inform their supervisor immediately.)

2.9.3 End Analysis and Finishing-Up

Upon completion of encapsulation, the operator must ensure the following is completed:

1. End Analysis' test, which is the same as the 'Start up Analysis' test described in section above.
2. QA officers must condition check the batch of capsules before the bags are sealed.
3. Collection of the 'Composite' sample for the Quality Control laboratory
4. Once QC approval has been obtained, weigh the batch and record the yield on the BMR. Send the completed capsule batch to the stores.

Note: Any remaining empty shells should be placed back into the supplier bag and sealed. The bag should then be placed back in the capsule box and again sealed. The Capsule shell part label is to be filled in with the weight of the remaining shells and returned to the warehouse. The weight of the returning shells is to be recorded in the BMR.

5. The BMR from the completed batch must be returned to the production office for checking by a MS.

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2.10 Capsule Filling and its Monitoring Process

You should give very special consideration to the choice of appropriate control procedures for the filling operations. Proper monitoring of the filling process is carried out through two specific methods. Firstly, it is carried out through process control and management of the encapsulations machine and related processual parameters which include the following:

- Standard fill weight
- Mechanical pellet filling chamber
- Machine Speed
- Capsule separation vacuum
- Compressed air pressure
- Mini-filled capsule sorter
- capsule Locking Length

Secondly, the monitoring ensures that the finalised, filled capsules meet the appropriate regulatory requirements which include weight variation, solubility, and disintegration which we shall discuss later in In-Process Control Tests for Capsules in subsection 2.11.

Now, let us further analyse some of the common defects found during the capsule production process (hard gelatin capsules) as illustrated below. In the following Table 2.2.

2.11 SOPs for In-Process Tests During


Fig 2.3. Visual Defects for 'Capsule Dependent Faults'


Fig 2.4. Other defects found in finished products

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Table 2.2 Common defects and causes in Capsule Production

PROBLEM	DESCRIPTION	CAUSE(S)
HOLE	An irregular opening in the cap or body	Supplier issue
SPLIT	A split in the film starting from the cap or body edge	Possible supplier issue Poor handling Misalignment of capsule segments Orienting issues
CRACKED	A cap or body with many splits	Possible supplier issue Poor handling Misalignment of capsule segments Orienting issues Closing issues
UNCUT BODY/CAP	An untrimmed cap or body	Supplier issue
FAILURE TO SEPARATE	Cap and body may not be separated properly	Possible supplier issue Blocked vacuum filter Insufficient vacuum. Blocked lower segment(s) High temperature during transport or storage
TRIMMING	A piece or the whole trimmed end of a cap or body inside a closed capsule	Possible supplier issue Poor handling Misalignment of capsule segments Orienting issues
DOUBLE DIP	Extra thick cap due to being dipped	Supplier issue
SHORT BODY	Cap length is 1mm less than the specified length	Supplier issue
JOINED IN LOCK	Capsule is in locked position	Supplier issue Orienting issues
DOUBLE CAP	A capsule with an additional cap covering the body end	Possible supplier issue Separation caused by excess vacuum during transport from capsule conveyor to capsule feed hopper. Separation during filling process, and capsule cap attaching to filled lower body during the subsequent polishing process
UNJOINED	A single cap or body	Possible supplier issue Separation caused by excess vacuum during transport from capsule conveyor to capsule feed hopper.

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THIN SPOT (CAP SHOULDER)	A thin area in the capsule shoulder which may rupture when the capsule is filled	Possible supplier issue Reduce vacuum at orienting station This issue may also result in 'cutting' of the end of the lower body.
MASHED	A mechanically damaged capsule that has been squashed flat	Poor handling Misalignment of capsule segments Orienting issues Use of scoop to remove excess filled product from bucket at the weighing stage
LONG CAP/BODY	Length of cap or body is 1mm more than the specified length	Supplier issue
PINCHED	The cap or body is damaged in the segments. Pinches >3mm are classified as a Major A defect.	Poor handling Misalignment of capsule segments Orienting issues
TELESCOPE	A closed capsule with a protruding piece of either cap or body produced by a double split.	Poor handling Misalignment of capsule segments Orienting issues Over-fill of powder/granulate in the filled capsule
LONG JOINED	Capsule not closed sufficiently to engage the pre-lock	Possible supplier issue such as capsule length Out-of- Specification (OOS) Incorrect setting at the Closing Station
DYE SPECK	A colored spot of pigment aggregate	Supplier issue
INVERTED END	A cap or body with the end pushed inwards	Possible supplier issue Poor handling Misalignment of capsule segments Orienting issues Problems at the discharge station
GREASE	Mould release aid spots on the inside of the capsule	Supplier issue
DAMAGED EDGE-LARGE	The edge of the cap is roughly trimmed. The imperfection at it is greatest >1mm into the specified length.	Supplier issue
SMALL PINCHED	Cap or body is damaged in the segment. Pinches <3mm are classified as a Major B defect.	Poor handling Misalignment of capsule segments Orienting issues
THIN SPOT	A thin area in the cap or body wall that may rupture when the capsule is filled.	Supplier issue

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Capsule Production

1. At all times when handling capsules, operators must take the opportunity to assess the capsule condition. If any defect is observed, highlight the issue immediately to the relevant Supervisor/QA officer.
2. Operators must take note of the testing frequency for the IPC tests outlined in the BMR as different markets have different testing requirements.
3. During encapsulation, there must be continuous condition checking of the buckets by QA officers, Manufacturing Supervisor or Team Leader. The encapsulation operator must not allow a build-up of unchecked buckets within the room, in case a capsule defect has gone unnoticed.
4. Ensure that the QA officer, Manufacturing Supervisor or Team Leader enters the capsule description onto the capsule 'Condition Check' page.
5. If for any reason the encapsulation equipment has to be cleaned during encapsulation, ensure to closely monitor the condition of the capsules produced after restarting the encapsulation machine.
6. During the encapsulation of the batch, the granulate/powder hopper volume must remain constant at all times until the IBC is almost empty. The batch must continue running until the granulate/powder in the powder/granulate feed hopper reaches the 'Level Control Initiator'. Never allow the encapsulation machine to run "dry" of granulate/powder.
7. Ensure that all steps are signed for in the BMR as each action takes place. Take care to ensure that IPC values entered into the BMR are within specification as outlined in the BMR.
8. If an Out of Specification (OOS) result is encountered during an 'In Process Check' (IPC), the machine must be stopped and the root cause must be established. The operator must use the following procedure:
 - Notify PS / QA Officer of the OOS result
 - Remove the ancillary equipment from the front of the machine.
 - Adjust the machine accordingly to achieve the previous encapsulation parameters.
 - Set up the ancillary equipment to the front of the machine when the encapsulation parameters are within the BMR specification.
 - If the previous encapsulation parameters cannot be achieved, then a QC officer or Production Supervisor must be notified and corrective actions must be taken before encapsulation can proceed. Any comments written in the BMR by the operator must be signed and dated.
9. Ensure that a composite sample is taken from each bucket compressed.
10. Ensure the samples are taken before the bucket is weighted and recorded in the Manufacturing BMR.
11. Ensure that each bucket is weighed to three decimal points with the weight entered into the BMR to three decimal points, and also on the corresponding bucket label.

You should note that in order to ensure the accuracy of this weighing operation, the floor balance is re-tared after every third bucket **if the product is not controlled. For controlled products, every bucket must be tared. If an additional composite sample is placed in the last bucket, the weight of the securitainer will be included in the tare weight.** Ensure that the balance used is set correctly and the level bubble for the balance is centred.

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12. Visually assess the top portion of each bucket for capsule defects after filling and prior to its removal from the floor balance.
13. Ensure that for those products that require desiccants that the correct quantity of desiccant bags are present and introduced to the tablet bucket according to the BMR.

2.12 Storage of Capsules

2.12.1 Empty capsules

Usually, empty capsules are shipped with a moisture contents between 13 to 16%. And you should ensure that this moisture contents is maintained. Also, avoid exposing the capsules to high temperatures or to cycles of high and low temperatures. There is also a tremendous volume of air inside the capsule that can extract or release moisture from the capsule. Maintain the area where you store the capsule at 15 to 25 ° C and 35 to 55% relative humidity (RH). Do not store empty capsules in freezers. The empty capsules are also very susceptible to damage because of the capsule walls are unsupported.

2.12.2 Filled capsules

The storage requirements for filled capsules are based on the product stability profile. The warehouse storage areas should be temperature-mapped and monitored to ensure the room temperature is under control and that product specific refrigeration/RH control requirements are maintained. Ideally, the storage areas would be equipped with alarms and there would be a product-specific product handling procedure that explains how to deal with out-of-limit temperature/humidity incidents. If the product is stored in bulk containers for a significant period of time prior to packaging, you should institute procedures to monitor bulk product stability.

2.12.3 General storage recommendations

Protect the capsules from direct sunlight by storing them away from windows and skylights. They should be stored away from radiators, heat registers, hot water pipes and team pipes.

Keep the capsules out from under potential sources of water condensation, such as water pipes, and keep pallet load off the floor.

Learning Activities

In groups, discuss how to avoid common problems of capsule production processes

Test Yourself

1. Which are the common encapsulation machines used in Jordan Pharmaceutical Manufacturing firms you visited?
2. Explain the IPCs employed during capsule production.

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